



AR App Development:

Pelvic Tumor Resection Guidance

Done by:

Alba Centeno, Noelia García, Antonio Sánchez and Marta Martín

TABLE OF CONTENTS

01

INTRODUCTION

02

OBJECTIVES

03

ROADMAP

04

RESULTS

05

CONCLUSION

06


FUTURE LINES

INTRODUCTION




BACKGROUND

Pelvic tumors pose significant challenges in surgical management due to their complex anatomical location and potential for aggressive growth.



INTRODUCTION



IMPORTANCE

Given the complexities associated with pelvic sarcomas and the challenges in achieving optimal resection margins while preserving limb function, tools such as AR technology can revolutionize surgical planning and execution.



OBJECTIVES

The objective of our project is to **develop an AR application** capable of calculating optimal margins and projecting the tumor along with its safety resection boundaries in real-time during surgery.

The aim is to provide the surgeon with a **visual and practical tool** that:

- ❑ Facilitates precise tumor resection
- ❑ Minimizes damage to surrounding structures
- ❑ Improves surgical and oncological outcomes for patients



SPECIFIC OBJECTIVES


Several key steps had to be considered for the accomplishment of the project.

This includes:

- ❑ Design and 3D printing of a **surgical guide** for attachment to the pelvic region.
- ❑ Definition of a **AR-based marker** for our application.
- ❑ Development of a **AR application for tumor visualization** on top of the pelvis.
- ❑ Added **functionalities** for clinical environment simulation.
 - Ability to customize for different patients and tumor types.
 - Option to show the optimal approach for resection.
 - Toggle options to show/hide the pelvic bone and tumor.
 - Option to display the coronal cross-section of the tumor.



ROADMAP



vuforia™



PELVIS MODEL

Model Utilized: Pelvis 3D printed Phantom

- Lent to us by the UC3M lab
- Guide attached to it
- Accurately replicate the anatomical structure of the human pelvis in the app

Virtual model:




Physical model:



SURGICAL GUIDE DESIGN

- Design **surgical guide** to attach to pelvis bone.
- Design **card holder** for AR marker.

Added side slot to introduce the card






 AUTODESK®
FUSION 360™

 AUTODESK
MESHMIXER




SURGICAL GUIDE DESIGN

- ❖ 3D print the combined model.
- ❖ Several tries
- ❖ Solution: changing printing parameters
 - Support structure to **Normal**
 - Build plate temperature to **60°**
 - Infill density to **10%**



MARKER DEFINITION

- ❖ Software: **vuforia**
- ❖ Requirement: universal marker
 - No personal identification
 - Easily identifiable by app
 - Available everywhere



vuforia™



PelviAR DEVELOPMENT

Main part of the project

Software: **Unity** & VSCode → 7 C# scripts



Design requirements:

1. Practical and **user-friendly** interface
2. Aid surgeons **prior** the operation
3. Visualize **resection margins** on patient




PelviAR FEATURES


- Dropdown for **patient selection**
- Slider for tumor **transparency**
- Button for superposition of real and **virtual model**



PelviAR FEATURES



- Button for **anatomical CT image**
- Calculation of **sphere boundary**
- **Input field** for personalized margin



LIMITATIONS

Lack of database




iOS implementation



LIMITATIONS

Margins projections



Tumor location



**Rudimentary margin
calculation**



FUTURE LINES

Solve limitations

- ❑ Accurate projection of margins and tumor
- ❑ Real 3D model from a CT.
- ❑ IOS implementation.

Advanced calculation margins algorithm

- ❑ Consider tumor size, severity and location
- ❑ Consider nearby structures (vascularization)



FUTURE LINES

Solve limitations

- ❑ Accurate projection of margins and tumor
- ❑ Real 3D model from a CT.
- ❑ iOS implementation.

Advanced calculation margins algorithm

- ❑ Consider tumor size, severity and location
- ❑ Consider nearby structures (vascularization)

Validation Method

- ❑ Tracking systems: evaluate the effectiveness of the app



CONCLUSION

FIRST STEP



APP

- Visualization
- Optimal margins calculation



IMPROVE PATIENT'S RECOVERY

Aids in the obtention of negative margins that have a positive impact on patient's recovery





THANKS

DO YOU HAVE ANY QUESTIONS?





PelviAR DEMO

